

Catalog Number AP 014

## Highlights:

- Test corn, cotton, canola or soybean
- High Sensitivity Protocol detects the PAT enzyme from the *pat* gene (PAT/*pat*) found in 0.5% T25 corn
- Rapid Protocol screens individual seeds or leaf samples for the presence of PAT/*pat*

## Contents of Kit:

- 1 antibody-coated 96-well plate
- PAT/*pat* Enzyme Conjugate
- 1 packet of Buffer Salts
- Substrate
- Stop Solution

## Intended Use

The QualiPlate Kit for LibertyLink PAT/*pat* is designed for the qualitative laboratory detection of phosphinothricin acetyl transferase enzyme (PAT) coded for by the *pat* gene in grain, leaf, or seed. Two assay protocols are presented. The High Sensitivity Protocol will detect the PAT enzyme from the *pat* gene (PAT/*pat*) found in 0.5% LibertyLink T25 corn, 40% Bt11(T14) corn, 20% Herculex™ I corn, 5% WideStrike™ cotton, 0.5% T45 LibertyLink canola, and 1.0% LibertyLink soybean, and requires 2.5 hours to run. The Rapid Protocol (one hour total) is intended for use in screening individual seeds or leaf samples for the presence of PAT/*pat*.

NOTE: LibertyLink PAT from the *bar* gene, at concentrations present in StarLink® corn, LL25 cotton, or InVigor™ canola, is not detected in either format.

## How the Test Works

This QualiPlate Kit is a “sandwich” Enzyme-Linked ImmunoSorbent Assay (ELISA). In the test, **corn, cotton, canola, or soybean** sample extracts are added to test wells coated with antibodies raised against PAT from the *pat* gene. Any residues present in the sample extract bind to the antibodies, and are then detected by addition of enzyme (horseradish peroxidase)-labeled PAT/*pat* antibody.

After a simple wash step, the results of the assay are visualized with a color development step; color development is proportional to PAT/*pat* concentration in the sample extract.

*Lighter color = Lower concentration*

*Darker color = Higher concentration*

## Materials not Provided

- distilled or deionized water for preparing Wash/Extraction Buffer
- glass bottles or flask plus graduated cylinder with 1 liter capacity for preparation and storage of Wash/Extraction Buffer
- test or centrifuge tubes for extraction of grain
- Tissue Extraction Kit: snap-cap tubes and pestles for extraction of leaf samples (EnviroLogix Cat. No. ACC 002, 100/package)
- centrifuge capable of 5000 x g
- disposable tip, adjustable air-displacement pipettes which will measure 50 and 100 microliters (µL)
- marking pen (indelible)
- tape or Parafilm®
- timer
- microtiter ELISA plate reader
- wash bottle, or microtiter plate or strip washer
- multi-channel pipette that will measure 50 and 100 µL
- racked dilution tubes for loading samples into the plate with a multi-channel pipette, or the equivalent
- orbital plate shaker (optional)



Prepare  
Wash/Extraction Buffer

### Preparation of Solutions

#### Wash/Extraction Buffer:

Add the contents of the packet of **Buffer Salts** to 1 liter of distilled or deionized water and stir to dissolve. Store refrigerated when not in use; warm to room temperature prior to assay. If more extraction buffer is needed, order item # P-3563 from Sigma Chemical Co. (St. Louis, MO), or prepare the equivalent.

## Choose the Assay Protocol

Choose the assay protocol (High Sensitivity or Rapid) that best suits the samples to be tested according to the table below. A “yes” in a column indicates that the indicated trait can be reliably detected in the assay protocol named. PAT/*pat* protein is expressed at levels too low to be reliably detected in the crops/tissues with a “no” statement.

Limits of detection (LOD) stated for bulk seed are in units of % by weight of PAT/*pat*-expressing seed in a bulk seed sample.

Crop	Event or Tradename	Leaf Rapid Protocol	Single Seed Rapid Protocol	Bulk ground seed High Sensitivity Protocol
Corn	LibertyLink T25	Yes	Yes	0.5% LOD
Corn	Herculex I	Yes	No	20% LOD
Corn	Bt11/T14	No	No	40% LOD
Cotton	WideStrike	No	No	5% LOD
Canola	LibertyLink T45	Yes	Yes	0.5% LOD
Soybean	LibertyLink	Yes	Yes	1.0% LOD

## Sample Preparation

Note: PAT/*pat* protein is not stable in solution. Extraction should be performed in 30 minutes or less. Do not extract overnight.

### High Sensitivity Protocol for Ground Grain/Seed

#### Corn:

This protocol requires that a small sample be analyzed. It is essential that this sample be well mixed and representative of the larger bulk. The test will detect 0.5% T25, 40% Bt11, or 20% Herculex I corn containing PAT/*pat*.

**Note:** Thorough mixing of the bulk grain sample and determination of an appropriate sampling plan are critical to the results of this testing, and are the responsibility of the user of this test kit. The USDA/GIPSA has prepared several guidance documents to address the issues involved in obtaining representative grain samples from static lots - such as trucks, barges, and railcars - and for taking samples from grain streams. Sampling plans should be chosen that best meet the needs of both the buyer and seller in terms of acceptable risks. Increasing the number of kernels in the sample and taking multiple samples will increase the likelihood of obtaining representative samples, and maximize the probability of detecting any contamination in the grain lot. For further information on USDA/GIPSA guidelines for obtaining representative samples and assessing detection probabilities for biotech grain, see the websites listed on at the left.

It is the responsibility of the user to ensure proper sampling and thorough mixing prior to analysis. Once representative samples have been obtained from the truck or container, they can be reduced in size using a splitter and uniformly ground and mixed.

**The finer the grind, the faster and more efficient the extraction.**

1. Weigh 5 grams of ground corn sample into a 30+ mL capacity vial or tube.
2. Add 12.5 mL of Wash/Extraction Buffer to each 5 gram sample. For all other grain sample sizes, add Wash/Extraction Buffer at the rate of 2.5 mL per gram of grain. Cap and shake vigorously by hand or vortex for 20-30 seconds. Let stand at room temperature for 10 to 30 minutes to extract.
3. The extracted samples/controls must be clarified by centrifuging the extract at 5000 x g for 5 minutes. Insert a pipette tip below any floating lipid layer and above the pellet to remove the clarified sample.

### Reference Websites:

[www.gipsa.usda.gov/fgis/handbook/gihbk1\\_inspec.aspx](http://www.gipsa.usda.gov/fgis/handbook/gihbk1_inspec.aspx)  
*USDA Grain Inspection Handbook, Book 1, Grain Sampling. This document provides a comprehensive overview of recommended sampling guidelines for static lots and grain streams. It reviews the various types of equipment and strategies that can be used to obtain a representative grain sample from different types of containers.*

[www.gipsa.usda.gov/fgis/biotech/sample2.htm](http://www.gipsa.usda.gov/fgis/biotech/sample2.htm)  
*Guidance document entitled Sampling for the Detection of Biotech Grains, which provides important statistical sampling considerations when testing for the presence of biotech grains. It covers the basis for making probability determinations in accepting lots based upon different assumptions with respect to sample size, number of samples, sample preparation, etc.*

[www.gipsa.usda.gov/fgis/biotech/sample1.htm](http://www.gipsa.usda.gov/fgis/biotech/sample1.htm)  
*Practical Application of Sampling for the Detection of Biotech Grains. This one-page application guide provides a table that gives sample sizes for selected lot concentrations and probability of rejecting the specified concentrations. It also provides a formula for making the calculation for other combinations.*

[www.gipsa.usda.gov/fgis/biotech/samplingplan1.xls](http://www.gipsa.usda.gov/fgis/biotech/samplingplan1.xls)  
*This website provides a simple to use Sample Planner (29k Excel Spreadsheet). The planner allows you to enter different assumptions in terms of sample size, number of samples, acceptable quality level and to determine the probability of accepting lots with given concentration levels. It also plots the probabilities in graph form for easy interpretation. Specific data can be saved for documentation and future analyses.*

### Canola:

1. Weigh out 2 grams of canola sample into a mesh bag (EnviroLogix Cat No. ACC 021). Use a rubber mallet to crush the seeds thoroughly. After crushing, run your fingers over the bag to make sure there are NO uncrushed seeds remaining.
2. Add 6.0 mL of Wash/Extraction Buffer to each 2 gram sample. Thoroughly massage the extraction buffer around in the bag until the entire sample is wet. Allow the sample to extract for 10 to 30 minutes at room temperature.
3. Remix the sample using your finger on the outside of the bag until you feel the sample is homogenous. Tilt the bag slightly and use your fingers to pull the extraction buffer over to one corner. Pipette out a sample (the sample will contain particles).
4. The extracted samples must be clarified by centrifuging the extract at 5000 x g for 5 minutes. Insert a pipette tip below any floating lipid layer and above the pellet to remove the clarified sample.

### Bulk Cottonseed:

1. Weigh out 5 grams ground cotton into 30+ mL capacity vial or tube.
2. Add 20 mL of Wash/Extraction Buffer to each 5 gram cotton sample. For all other cotton sample sizes, add Wash/Extraction Buffer at the rate of 4.0 mL per gram of cotton. Cap and shake vigorously by hand or vortex for 20-30 seconds. Let stand at room temperature for 10 to 30 minutes to extract.
3. The extracted samples must be clarified by centrifuging the extract at 5000 x g for 5 minutes. Insert a pipette tip below any floating lipid layer and above the pellet to remove the clarified sample.

### Soybean Grain:

1. Weigh out 5 grams ground soybean into 30+ mL capacity vial or tube.
2. Add 25 mL of Wash/Extraction Buffer to each 5 gram soy sample. For all other soy sample sizes, add Wash/Extraction Buffer at the rate of 5.0 mL per gram of soy. Cap and shake vigorously by hand or vortex for 20-30 seconds. Let stand at room temperature for 10 to 30 minutes to extract.
3. The extracted samples must be clarified by centrifuging the extract at 5000 x g for 5 minutes. Insert a pipette tip below any floating lipid layer and above the pellet to remove the clarified sample.

### Rapid Protocol for Screening Single Seeds or Leaf Punches:

#### Individual seeds:

1. Crush seeds: Corn kernels, cotton seeds, or single soybeans may be placed in a plastic bag and smashed with a hammer, then transferred to a tube; or a seed crusher/48-well plate combination may be used (Hypure #HSC-100, PerkinElmer, with Costar plate #3548, Corning Life Sciences, or equivalent). Check to be sure that all seeds have been crushed. Take extreme care not to cross-contaminate between seed samples. Canola seeds can be crushed by placing them into the wells of an uncoated 96-well plate and mashing each seed with the flat end of the pestle (EnviroLogix Cat No. ACC 002). Wipe the pestle clean on a damp paper towel between each seed.
2. Add 0.5 mL of Wash/Extraction Buffer to each crushed corn seed, 1 mL to each crushed cotton seed, 0.2 mL to each crushed canola seed, or 1 mL to each crushed soybean. Mix for at least 30 seconds, then allow particles to settle.



Allow all reagents to reach room temperature before beginning



Leaf punch



Extract sample



Centrifuge to clarify sample extract  
(High Sensitivity protocol only)

### Leaf samples:

1. Take a single leaf punch of approximately 10 mm diameter, using a micro-tube cap. Mash the leaf tissue with a pestle matched to the micro-tube (EnviroLogix Tissue Extraction Kit).
2. Add 0.5 mL of Wash/Extraction Buffer per leaf punch. Mix for at least 30 seconds, then allow particles to settle. Take extreme care not to cross-contaminate between leaf samples.

Note: Leaf samples must be run in the assay on the same day they are punched.

## How to Run the Assay

- Read all of these instructions before running the kit.
- Allow all reagents to reach room temperature before beginning (at least 30 minutes with un-boxed plates and reagents at room temperature - do not remove plate from bag with desiccant until it has warmed up).
- Organize all sample extracts, and pipettes so that Step 1 can be performed in 15 minutes or less. The use of a multi-channel pipette is strongly recommended.
- Use a disposable-tip air-displacement pipette and a clean pipette tip to add each Calibrator and sample extract to the wells. Conjugate, Substrate, and Stop Solution may be added in the same manner; alternatively, use a repeating pipette with a disposable tip on the end of the Combitip for these three reagents.
- Use the well identification markings on the plate frame to guide you when adding the samples and reagents. For this qualitative assay, duplicate wells of the Wash/Extraction Buffer Blank (BL), user-supplied known negative sample (Negative Control, NC), and user-supplied known positive sample (Positive Control, PC), along with 90 sample extracts (S) in single wells may be run on one plate. (See the Qualitative Assay Example Plate Layout - Figure 1A, page 6).


1. Add **50 µL** of **Enzyme Conjugate** to each well, followed immediately by **50 µL** of **Wash/Extraction Buffer Blank (BL)**, **50 µL** of clarified **Positive** and **Negative Control sample extracts (PC, NC)**, and **50 µL** of each clarified **sample extract (S)** to their respective wells, as shown in the Example Plate Layout (Figure 1A).

**NOTE:** In order to minimize setup time it is recommended that a multi-channel pipette be used in steps 1, 5, and 7.


2. Thoroughly mix the contents of the wells by moving the plate in a rapid circular motion on the benchtop for a full 20-30 seconds. Be careful not to spill the contents!
3. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature according to the table below. If an orbital plate shaker is available, shake plate at 200 rpm.

<b>HIGH SENSITIVITY PROTOCOL</b>	<b>2 hours</b>
<b>RAPID PROTOCOL</b>	<b>45 minutes</b>

4. After incubation, carefully remove the covering and vigorously shake the contents of the wells into a sink or other suitable container. Flood the wells completely with Wash/Extraction Buffer, then shake to empty. Repeat this wash step three times. Alternatively, perform these four washes (300 µL/well) with a microtiter plate or strip washer. Slap the inverted plate on a paper towel to remove as much liquid as possible.



*Add Conjugate, controls and samples*




*Mix plate*



*Incubate*



*Wash plate*



*Read plates in a Plate Reader within 30 minutes of the addition of Stop Solution*

5. Add **100 µL** of **Substrate** to each well.
6. Thoroughly mix the contents of the wells, as in step 2. Cover the wells with new tape or Parafilm and incubate at ambient temperature according to the table below. Use orbital shaker if available.

<b>HIGH SENSITIVITY PROTOCOL</b>	<b>30 minutes</b>
<b>RAPID PROTOCOL</b>	<b>15 minutes</b>

**Caution: Stop Solution is 1.0N Hydrochloric acid. Handle carefully.**

7. Add **100 µL** of **Stop Solution** to each well and mix thoroughly. This will turn the well contents yellow.

**NOTE:** Read the plate within 30 minutes of the addition of Stop Solution.

## How to Interpret the Results

### Spectrophotometric Measurement

1. Set the wavelength of the microtiter plate reader to 450 nanometers (nm). (If it has dual wavelength capability, use 600, 630 or 650 nm as the reference wavelength.)
2. Set the plate reader to blank on the Wash/Extraction Buffer Blank wells (this should automatically subtract the mean optical density (OD) of the Blank wells from each control and sample OD). If the reader cannot do this, it must be done manually.

*Note:* The mean OD of the BLANK wells should not exceed 0.2.

### Interpret the Results

A sample is generally considered to be positive if the blank-subtracted OD is greater than or equal to 0.2. For best results, a known negative sample as well as a low-level positive sample should be run as controls in each assay to help judge the results (a conventional and a 0.5% T25 ground corn sample, for example).

**Figure 1A. Example of a typical Qualitative assay setup.**

	1	2	3	4	5	6	7	8	9	10	11	12
A	BL	S6	S14	S22	S30	S38	S46	S54	S62	S70	S78	S86
B	NC	S7	S15	S23	S31	S39	S47	S55	S63	S71	S79	S87
C	PC	S8	S16	S24	S32	S40	S48	S56	S64	S72	S80	S88
D	S1	S9	S17	S25	S33	S41	S49	S57	S65	S73	S81	S89
E	S2	S10	S18	S26	S34	S42	S50	S58	S66	S74	S82	S90
F	S3	S11	S19	S27	S35	S43	S51	S59	S67	S75	S83	BL
G	S4	S12	S20	S28	S36	S44	S52	S60	S68	S76	S84	NC
H	S5	S13	S21	S29	S37	S45	S53	S61	S69	S77	S85	PC

## Precautions and Notes

- Store all Plate Kit components at 4°C to 8°C (39°F to 46°F) when not in use.
- Do not expose Plate Kit components to temperatures greater than 37°C (99°F) or less than 2°C (36°F).
- Allow all reagents to reach ambient temperature (18°C to 27°C or 64°F to 81°F) before use.



- Do not use kit components after the expiration date.
- Do not use reagents or plates from one Kit with reagents or plates from a different Kit.
- Do not expose Substrate to sunlight during pipetting or while incubating in test wells.
- The assay has been optimized to be used with the protocol provided in the kit. Deviation from this protocol may invalidate the results of the test.
- As with all tests, it is recommended that results be confirmed by an alternate method when necessary.
- Observe any applicable regulations when disposing of samples and kit reagents.
- Use extreme caution to prevent sample-to-sample cross-contamination with samples, fluids, or disposables.



**For Technical Support  
Contact Us At:**

**EnviroLogix**  
500 Riverside Industrial  
Parkway  
Portland, ME 04103-1486  
USA

**Tel: (207) 797-0300**  
**Toll Free: 866-408-4597**  
**Fax: (207) 797-7533**

e-mail:  
[info@envirologix.com](mailto:info@envirologix.com)

website:  
[www.envirologix.com](http://www.envirologix.com)



## LIMITED WARRANTY

EnviroLogix Inc. (“EnviroLogix”) warrants the products sold hereunder (“the Products”) against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product’s printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix’ sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix’ employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

*Parafilm is a registered trademark of American Can Corporation*  
*LibertyLink is a trademark of Bayer*  
*StarLink is a trademark of Aventis CropScience*  
*Herculex and WideStrike are trademarks of Dow AgroSciences LLC*  
*EnviroLogix, the EnviroLogix logo, and QualiPlate are trademarks of EnviroLogix Inc.*

© EnviroLogix 2015



**Safety Data Sheet**  
According to OSHA 29CFR 1910.1200

**SECTION 1. Identification of the substance/mixture and of the company/undertaking**

**1.1 Product identifier**  
Trade name: **Wash Buffer Salts**  
Part number: 50-0091, 10099

**1.2** Relevant identified uses of the substance or mixture and uses advised against application of the substance or the preparation: Laboratory chemicals

**1.3 Details of the supplier of the safety data sheet**  
Manufacturer/Supplier: EnviroLogix Inc., 500 Riverside Industrial Pkwy, Portland ME 04103, USA (207) 797-0300 (207) 797-0300 Technical Service

**1.4 Emergency telephone number:**

**SECTION 2. Hazards identification**

**2.1 Classification of the Substance or Mixture:**  
Classification according to OSHA 29CFR 1910.1200 (Hazard Communication): Not a hazardous substance or mixture

**2.2 Label Elements:**  
None required according to 29CFR 1910.1200

Other indications: None

**2.3 Additional information:**  
No other information

**SECTION 3. Composition/information on ingredients**

**3.1 Mixture:** Potassium acid phosphate  
Synonyms: PBS

Hazardous Components	Chemical name	CAS No	EC No	Amount (%)	Classification
	Potassium Chloride	7447-40-7	231-211-8	1-5 %	Aquatic Acute 3; Aquatic Chronic 3, H412

Based on the amount of hazardous ingredients in this product, it is not considered hazardous according to 29CFR 1910.1200

**SECTION 4. First aid measures**

**4.1 Description of first aid measures:**  
After inhalation: Supply fresh air, consult doctor in case of breathing difficulties.  
After skin contact: Flush skin with plenty of water for at least 15 minutes. Remove contaminated clothing. Seek medical attention if irritation develops.  
After eye contact: Rinse opened eye for several minutes under running water. Seek medical attention if irritation develops.  
After swallowing: If swallowed, consult with medical staff or poison control center to determine if any immediate response or follow up actions are recommended. Never give anything by mouth to an unconscious person.

**4.2 Most important symptoms and effects, both acute and delayed:** None

**4.3 Indication of any immediate medical attention and special treatment needed:** No special treatment is required

**SECTION 5. Firefighting measures**

**5.1 Extinguishing media:**  
Suitable extinguishing agents: CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture:**  
Carbon oxides, Oxides of Phosphorus, Potassium, Sodium, Hydrogen Chloride gas

**5.2 Advice for firefighters:**  
Wear protective equipment appropriate for fire conditions including respiratory protective gear

**SECTION 6. Accidental release measures**

**6.1 Personal precautions, protective equipment and emergency procedures:**  
Use PPE, avoid dust formation, ensure adequate ventilation, avoid breathing dust

**6.2 Environmental precautions:**  
Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge to the environment must be avoided.

**6.3 Methods and material for containment and clean up:**  
Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable closed containers for disposal

**6.4 Reference to other sections:**  
For safe handling refer to Section 7; For information on PPE refer to Section 8. For disposal, refer to Section 13.

**SECTION 7. Handling and storage**

**7.1 Precautions for safe handling:**  
Practice good chemical hygiene when handling. Avoid contact with eyes, skin and clothing. Prevent formation of dust.

**7.2 Conditions for safe storage, including any incompatibilities:**  
Keep containers closed, store in a dry, well ventilated space.

**7.3 Specific end use(s):**  
Apart from the uses mentioned in section 1.2, no other end uses are stipulated.

**SECTION 8. Exposure controls/personal protection**

**8.1 Control parameters:**  
Components with workplace control Parameters: Contains no substances with occupational exposure limit values

**8.2 Exposure controls:**  
8.2.1 Appropriate engineering controls: Ensure eyewash and safety shower are nearby; provide ventilation if necessary

8.2.2 Personal Protective Equipment:  
Eyes: Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection equipment are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals

Hands: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Respiratory protection: Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Body: Use body protection relative to its type and amount of material being handled

8.2.3 Environmental controls: Sweep or wipe up spills, do not allow into sewers or drains

**SECTION 9. Physical and chemical properties**

**9.1 Information on basic physical and chemical properties:**

a) Appearance:	White powder.
b) Odor:	None
c) Odor Threshold:	No data available
d) pH:	7-8
e) Melting point/freezing point:	No data available
f) Boiling point/boiling range:	No data available
g) Flash point:	No data available
h) Evaporation rate:	No data available
i) Flammability (solid, gaseous):	No data available
j) Upper/lower flammability or explosive limits:	No data available
k) Vapor pressure:	No data available
l) Vapor density:	No data available
m) Relative density:	No data available
n) Solubility(ies):	Water soluble
o) Partition coefficient: n-Octanol/water:	No data available
p) Auto-ignition temperature:	No data available
q) Decomposition temperature:	No data available
r) Viscosity:	No data available
s) Explosive properties:	No data available
t) Oxidizing properties:	No data available

**9.2 Other information:** No further relevant information available.

**SECTION 10. Stability and reactivity**

**10.1 Reactivity:** No data available

**10.2 Chemical stability:** Stable under normal recommended storage conditions.

**10.3 Possibility of hazardous reactions:** No data available

**10.4 Conditions to avoid:** No data available

**10.5 Incompatible materials:** Strong oxidizing agents and strong acids.

**10.6 Hazardous decomposition products:** No data available

**SECTION 11. Toxicological information**

Acute toxicity: No data available

Inhalation: No data available

Dermal: No data available

Skin corrosion/irritation: No data available

Serious eye damage: No data available

Respiratory or skin sensitization: No data available

Mutagenicity and toxicity for reproduction: No data available

Carcinogenicity: No component of this product at levels greater than 0.1 % is identified as probable, possible, or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA.

**SECTION 12. Ecological information**

**12.1 Toxicity:** No data available

**12.2 Persistence and degradability:** No data available

**12.3 Bio accumulative potential:** No data available

**12.4 Mobility in soil:** No data available

**12.5 Results of PBT and vPvB assessment:** Not available as a chemical safety assessment, not required/not conducted.

**12.6 Other adverse effects:** No data available

**SECTION 13. Disposal considerations**

Dispose of excess or unused product in accordance with Local, State and Federal regulations. Contact a licensed professional waste disposal service to dispose of this material.

**SECTION 14. Transport information**

**14.1 UN-number** (DOT, ADR, ADN, IMDG, IATA): Not dangerous goods

**14.2 UN proper shipping name** (DOT, ADR, ADN, IMDG, IATA): Not dangerous goods

**14.3 Transport hazard classes** (DOT, ADR, ADN, IMDG, IATA): Not applicable

**14.4 Packing group** (DOT, ADR, IMDG, IATA): Not applicable

**14.5 Environmental hazards:** Not applicable

**14.6 Special precautions for user:** Not applicable

**14.7 Transport in bulk according to Annex II of MARPOL 73/78:** Not applicable

**SECTION 15. Regulatory information**

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

**US Federal Regulations**  
SARA Section 302 (Extremely Hazardous Substances): Not listed  
Clean Air Act: Not listed  
Clean Water Act: Not listed  
OSHA: Not listed

**US State Regulations**  
Massachusetts Right to Know: Disodium Hydrogenorthophosphate CAS No 7558-79-4 Rev Date: 2007-03-01  
California Prop. 65 Components: Contains no chemicals known to cause cancer, birth defects, or reproductive harm.

**15.2 Chemical Safety Assessment:** Not carried out

**SECTION 16. Other information**

**Hazard Code**  
H412 Harmful to aquatic life with long lasting effects

*This information is true based on our present knowledge. However, EnviroLogix makes no representation of its accuracy or completeness. Persons receiving this information must exercise their independent judgment in determining the product's safety and suitability for its intended use. This document shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

EFES Department  
EnviroLogix Inc.





**Material Safety Data Sheet**  
OSHA 29CFR 1910.1200

**SECTION 1. Identification of the substance/mixture and of the company/undertaking**

<b>1.1 Product identifier</b>	<b>Stop Solution</b>
Trade name:	L.O.N HCl
Synonyms:	10825, 10827, 10828, 11193, 11776 (XGD007)
Part number:	Laboratory chemicals
<b>1.2 Relevant identified uses of the substance or mixture and uses advised against application of the substance / the preparation :</b>	
<b>1.3 Details of the supplier of the safety data sheet</b>	Envirol ogix Inc., 500 Riverside Industrial Pkwy, Portland ME, 04103, USA Phone: (207) 7974300
<b>1.4 Emergency telephone number:</b>	(207) 797-0300 Technical Service

**SECTION 2. Hazards identification**

<b>2.1 Classification of the substance or mixture</b>	<b>Hazard Classes</b>
Classification according to OSHA 29 CFR 1910.1200	Metal Corrosive (Cat. 1) H290 Skin Irritation (Cat 2) H315 Serious Eye damage (Cat. 1) H318
<b>2.2 Label elements</b>	
Labeling according to OSHA 29CFR 1910.1200	
Hazard pictograms :	
Signal word :	Warning
Hazard statements:	H290 May be corrosive to metals H315 Causes skin irritation H318 Causes serious eye damage
Precautionary statements:	P281 Use personal protective equipment as required P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.
<b>2.3 Other Statements</b>	None

<b>6.3 Methods and material for containment and cleanup:</b>	Absorb in paper towel and discard in appropriate waste. Clean with water afterwards. Large spills may be neutralized with dilute solutions of sodium carbonate or calcium oxide.
<b>6.4 References to other sections:</b>	For safe handling refer to Section 7. For information on PPE refer to Section 8. For disposal refer to Section 13.

**SECTION 7. Handling and storage**

<b>7.1 Precautions for safe handling:</b>	Practice good chemical hygiene when handling. Avoid contact with eyes, skin, and clothing.
<b>7.2 Conditions for safe storage, including any incompatibilities:</b>	Store in tightly closed, non-metal container, in a corrosive compatible area. Prevent direct sunlight and heat. Store in well aired storage rooms.
<b>7.3 Specific end use(s):</b>	Apart from the uses mentioned in section 1.2., no other specific uses are stipulated.

**SECTION 8. Exposure controls/personal protection**

<b>8.1 Exposure limits:</b>	Components with limit values that require monitoring at the workplace:									
	<table border="1"> <thead> <tr> <th>Hydrogen Chloride</th> <th>European (Commission directive 96/94)</th> <th>USA (OSHA)</th> </tr> </thead> <tbody> <tr> <td></td> <td>8hr TWA = 5 ppm (7.5 mg/m<sup>3</sup>)</td> <td>Ceiling Limit = 5 ppm (7.5 mg/m<sup>3</sup>)</td> </tr> <tr> <td></td> <td>STEL = 10 ppm (15 mg/m<sup>3</sup>)</td> <td></td> </tr> </tbody> </table>	Hydrogen Chloride	European (Commission directive 96/94)	USA (OSHA)		8hr TWA = 5 ppm (7.5 mg/m <sup>3</sup> )	Ceiling Limit = 5 ppm (7.5 mg/m <sup>3</sup> )		STEL = 10 ppm (15 mg/m <sup>3</sup> )	
Hydrogen Chloride	European (Commission directive 96/94)	USA (OSHA)								
	8hr TWA = 5 ppm (7.5 mg/m <sup>3</sup> )	Ceiling Limit = 5 ppm (7.5 mg/m <sup>3</sup> )								
	STEL = 10 ppm (15 mg/m <sup>3</sup> )									
<b>8.2 Exposure Controls:</b>	Facilities using this mixture should be equipped with an eyewash and safety shower. Use general or local exhaust ventilation to keep airborne concentrations below permissible exposure limits.									
8.2.1 Engineering controls										
8.2.2 General protective and hygienic measures:	The usual precautionary measures should be adhered to when handling chemicals.									
Eye Protection:	Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection regulations are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals.									
Hand Protection:	Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.									
Breathing Equipment:	Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).									
8.2.3 Environmental exposure controls:	Contain spills, do not allow into environment									

**SECTION 3. Composition/information on ingredients**

<b>3.2 Mixture</b>	Aqueous solution 1N Hydrochloric Acid (1N HCl, 3% HCl)			
Chemical name	Amount (%)	CAS No		Classification According to OSHA 29CFR 1910.1200
		EC No		
Hydrochloric acid	1-4 %	7647-01-0		Hazard Classification May be Corrosive to Metals Causes Skin Irritation Causes Serious Eye Damage
		231-595-7		

**SECTION 4. First aid measures**

<b>4.1 Description of first aid measures</b>	
After inhalation :	In case of inhalation: Remove to fresh air. If not breathing give artificial respiration. Get medical attention immediately.
After skin contact :	In case of skin contact: Remove contaminated clothing and shoes immediately. Wash affected area with mild soap or detergent for at least 10 minutes or until no evidence of chemical remains.
After eye contact :	In case of eye contact, immediately flush eyes with plenty of water for at least 15 minutes. Lifting eyelids occasionally, until no evidence of chemical remains. Get medical attention immediately.
After swallowing :	In case of ingestion, DO NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Call a physician immediately.
<b>4.2 Most important symptoms and effects, both acute and delayed:</b>	May cause skin irritation and eye damage
<b>4.3 Indication of any immediate medical attention and special treatment needed:</b>	DO NOT use sodium bicarbonate in an attempt to neutralize the acid.

**SECTION 5. Firefighting measures**

<b>5.1 Extinguishing media:</b>	CO <sub>2</sub> , extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
<b>5.2 Special hazards arising from the substance or mixture:</b>	Hydrogen Chloride gas
<b>5.3 Advice for firefighters:</b>	Wear protective gear appropriate for fire conditions including respiratory protective gear.

**SECTION 6. Accidental release measures**

<b>6.1 Personal precautions, protective equipment and emergency procedures:</b>	In the case of spilled mixture wear gloves to prevent skin contact. In the case of a large spill, additional protection is recommended.
<b>6.2 Environmental precautions:</b>	Do not discharge mixture to sewer system or waterways.

**SECTION 9. Physical and chemical properties**

<b>9.1 Information on basic physical and chemical properties:</b>	
a) Appearance:	Clear liquid, colorless to slight yellow.
b) Odor:	Pungent (slight)
c) Color Threshold:	No Data Available
d) pH:	pH 1
e) Melting point/freezing point:	No Data Available
f) Boiling point/Boiling range:	No Data Available.
g) Flash point:	Not applicable.
h) Evaporation rate:	0.36 (Water) compared with n-Butyl Acetate = 1
i) Flammability (solid, gaseous):	No Data Available
j) Upper/lower flammability or explosive limits:	No Data Available
k) Vapor pressure:	No Data Available
l) Vapor density:	No Data Available
m) Relative density:	No Data Available
n) Solubility(ies):	Fully miscible, water.
o) Partition Coefficient: n-Octanol/water:	No Data Available
p) Auto-ignition temperature:	No Data Available
q) Decomposition temperature:	No Data Available
r) Viscosity:	No Data Available but should be similar to that of water
s) Explosive properties:	No Data Available.
t) Oxidizing properties:	No Data Available
<b>9.2 Other information:</b>	No further relevant information available.

**SECTION 10. Stability and reactivity**

<b>10.1 Reactivity:</b>	No data available
<b>10.2 Chemical Stability:</b>	Stable under normal temperatures and pressures.
<b>10.3 Possibility of hazardous reactions:</b>	Under normal conditions of storage and use, hazardous reactions will not occur.
<b>10.4 Conditions to avoid:</b>	No specific data
<b>10.5 Incompatible materials:</b>	Metals, Alkali metals, bases, Amines.
<b>10.6 Hazardous decomposition products:</b>	Under normal conditions of storage and use, hazardous decomposition products should not be produced.

**SECTION 11. Toxicological information**

<b>Information on Toxicological Effects</b>													
Acute effects (toxicity tests):	<table border="1"> <thead> <tr> <th>7647-01-0 HCl</th> <th>Effect Dose</th> <th>Species</th> </tr> </thead> <tbody> <tr> <td>Acute oral toxicity</td> <td>LD50=900mg/kg</td> <td>rabbit</td> </tr> <tr> <td>Acute dermal toxicity</td> <td>No data</td> <td></td> </tr> <tr> <td>Acute inhalative toxicity</td> <td>LC50 = 3124 mg/L</td> <td>rat</td> </tr> </tbody> </table>	7647-01-0 HCl	Effect Dose	Species	Acute oral toxicity	LD50=900mg/kg	rabbit	Acute dermal toxicity	No data		Acute inhalative toxicity	LC50 = 3124 mg/L	rat
7647-01-0 HCl	Effect Dose	Species											
Acute oral toxicity	LD50=900mg/kg	rabbit											
Acute dermal toxicity	No data												
Acute inhalative toxicity	LC50 = 3124 mg/L	rat											
Sensitization:	No sensitizing effects known												
CMR (carcinogenicity, mutagenicity and toxicity for reproduction) effects:	No CMR effects.												
Additional toxicological information:	No Additional Information												

**SECTION 12. Ecological information**

<b>12.1 Toxicity:</b>			
Aquatic toxicity (1N HCl)	Effect dose	Exposure time	Species
Acute fish toxicity	LC50=826 mg/L	96h	Lemiscus idus
Acute daphnia toxicity	No data		
Acute algae toxicity	No data		