PREPARO DA AMOSTRA

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Prepare the sample

Collect and grind a representative corn sample such that 95% passes through a 20-mesh sieve - CORRECT

Example of a sample ground too coarse which will result in improper extraction – INCORRECT

Example of a sample ground too fine which will impede the flow of extract up the strip – INCORRECT

Test procedure

1. Add a sub-sample of 25g to the cup


2. Add one pouch of EB17, then 75 mL of water (distilled, deionized or bottled)

3. Shake immediately for 10 seconds to dissolve the pouch, then for 1 minute on a mechanical shaker or 2 minutes by hand

4. Filter the sample into a clean vessel with an approved coffee filter; move the filter to the side to access filtered extract

5. Pipette 100 μ L of DB2 into the reaction vial, then 100 μ L of filtered extract. Pipette up and down to homogenize well.

6. Place the QuickTox strip in the vial, wait 4 minutes. Remove strip immediately and interpret.

Test Interpretation

The appearance of the Test Line (regardless of its intensity) indicates that the sample contains less than 20 ppb, thus negative. If no Test Line develops within 4 minutes, it indicates that the sample contains more than 20 ppb, thus, is positive.