# QuickStix™ Kit for QuickScan Cry1Ac Soy Bulk Soybeans

Catalog Number AQ 433 BGBR

Part #11572

**Highlights:**

- Results in 10 minutes or less
- Available as 100-strip individual kits or bulk-packaging

**Contents of Kit:**

- 100 QuickStix Strips packed in moisture-resistant canisters
- 100 transfer pipettes
- 100 sample cups

**Items Not Provided:**

- Blender for sample prep: Oster® Sunbeam blender model #4000 series
- Graduated cylinder
- Tap water
- QuickScan System

Grind sample

**Soybean Extraction**

Grams of beans x 3.5 = mL tap water

For example: (1000 x 0.15 average seed weight) = 150g

150g x 3.5 = 525 mL tap water

**Intended Use**

The EnviroLogix QuickStix Kit for QuickScan - Cry1Ac Bulk Soybeans is designed to detect and quantify the presence of Cry1Ac protein at the levels typically expressed in MON 87701 soybean. The detection level of these QuickStix Strips varies depending on the expression level of the soybean; this equates to:

High expressor	0.1%	1 bean in 1000
Medium expressor	0.15%	1 bean in 650
Low expressor	0.3%	1 bean in 333

The test may be used qualitatively without the use of QuickScan. A visible pink test line will appear if the sample contains 0.3% Cry1Ac or more.

**How the Test Works**

In order to detect the Cry1Ac protein expressed in the MON 87701 soybeans, the sample must first be extracted to solubilize the protein. Each QuickStix Strip has an absorbent pad at each end. The protective tape with the arrow indicates the end of the strip to insert into the reaction cup. The sample will travel up the membrane strip and be absorbed into the larger pad at the top of the strip. The portion of the strip between the protective tape and the absorbent pad at the top of the strip is used to view the reactions as described under “Interpreting the Results”. Please avoid bending the strips.

**Sample Preparation**

**Step 1: Determine Size of Sub-samples**

1. Determine the average individual weight of seeds/grain being tested, weighing 100 grains and dividing by 100.
2. Calculate the total weight of the seeds/grains to be tested by multiplying the number of seeds/grains x average seed/grain weight. A sample size of 1000 seeds/grain is recommended; other sample sizes may also be used.

**Step 2: Prepare the sample extract**

3. Weigh or count beans into the Oster blender (or appropriate vessel).
4. Grind sample on high speed for 30 seconds.
5. Add the volume of buffer calculated by the formula (left).
6. Pulse blender three times to mix well.
7. Let sample settle for 20-30 seconds then draw off liquid at that time. Avoid pulling up particles with the transfer pipette.
8. Transfer approximately 2 mL of the liquid extract from above the settled soybean solids into the auto analyzer cup. The level of extract should be to the ridge of the cup.
9. To prevent cross-contamination, thoroughly clean blender parts and jars of dust and residue prior to preparation of a second sample. Use a new transfer pipette and cup for each sample.

Note: For instructions on running this test along with the QuickStix Kit for QuickScan-Roundup Ready bulk soybeans, see instructions under Precautions and Notes.

**How to Run the QuickStix Strip Test**

1. Allow refrigerated canisters to come to room temperature before opening. Remove the QuickStix Strips to be used. Avoid bending the strips. Reseal the canister immediately.

Fill cup to ridge with extract

Add QuickStix Strip,  
wait 10 minutes

Any pink Test Line  
is considered positive

Remove arrow tape, place  
QuickStix in QuickScan Reader

- Place the strip into the sample cup containing 2 mL of the liquid soybean extract. The sample will travel up the strip.
- Allow the strip to develop for 10 minutes before making final assay interpretations. Positive sample results may become obvious much more quickly.
- Remove strip immediately; cut off and discard the bottom section of the strip covered by the arrow tape; and place QuickStix Strip(s) in the QuickScan Reader – strips must be read immediately after cutting, while still wet. Alternatively, strips may be read visually at this point. Do not save strips.

**NOTE:** Use extreme caution to prevent sample-to-sample cross-contamination with grain/seeds, fluids, materials, or disposables.


### Interpreting the Results

Development of the Control Line within 10 minutes indicates that the strip has functioned properly. Any strip that does not develop a Control Line should be discarded, and the sample re-tested using another strip.

Results are scanned and interpreted quantitatively with the QuickScan System. Place QuickStix into the carrier, slide in, and press “Read Test” on the screen. QuickScan will return a result as “% GMO” or “<LOD” (less than the Limit of Detection). Please consult the QuickScan User Manual for details.

**NOTE:** Because of the variability in expression, it is possible to have a very dark Test Line and test area that can cause QuickScan to return a result “E10” – this should be considered a 100% positive result.

The test may be used qualitatively without the use of QuickScan. A visible pink test line will appear if the sample contains 0.3% Cry1Ac or more.

### Kit Storage

QuickStix can be stored at room temperature, or refrigerated for a longer shelf life. Note the shelf life on the kit box for each storage temperature. The kit may be used in field applications; however, prolonged exposure to high temperatures may adversely affect the test results. Do not open the desiccated canister until ready to use the test strips.

### Precautions and Notes

- This kit is designed to give a quantitative result using the QuickScan System and is not intended to be visually interpreted.
- This product is currently not applicable for use in any other crop.
- This product will not perform as expected if the seed/grain is treated.
- As with all tests, it is recommended that results be confirmed by an alternate method when necessary.
- The assay has been optimized to be used with the protocol provided in the kit. Deviation from this protocol may invalidate the results of the test.
- The results generated through the proper use of this diagnostic tool reflect the condition of the working sample directly tested. Extrapolation as to the condition of the originating lot, from which the working sample was derived, should be based on sound sampling procedures and statistical calculations which address random sampling effects, non-random seed lot sampling effects and assay system uncertainty. A negative result obtained when properly testing the working sample does not necessarily mean the originating lot is entirely negative for the analyte or protein in question.
- A negative result with this kit does not mean that the sampled tissue has not been otherwise genetically modified.
- Protect all components from hot or cold extremes of temperature when not in use.

- CAUTION: Tightly closed containers of soy extract, if left sitting for several hours, may ferment and cause the lid or container to burst. Dispose of extract when testing is complete.

**Running this Kit Along With QuickStix Kit for QuickScan-Roundup Ready**

Testing vials may vary – be sure to use the vials supplied with the kit to run the strips in that kit.

A. Weigh or count beans into Oster blender, grind		
B. Add 3.5 mL/g of water, pulse blender 3 times to mix		
C. Let sample settle 20-30 seconds, then transfer 2 mL to Cry1 auto analyzer cup		
D. Add Cry1 strip and start timer	Test for 10 minutes	H. After testing time, remove each strip, cut off lower portion and discard, and place test portion into carrier; slide in and press Read.
E. Add another 3.5 mL/g of water, pulse blender to mix again		
F. Allow sample to settle 20-30 seconds then transfer extract to RR auto analyzer cup per the Product Insert		
G. Add RR strip and start timer	Test for 5 minutes	

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AQ 433 BGBR

Usar tubo de 2 mL

EnviroLogix Inc. (“EnviroLogix”) warrants the products sold hereunder (“the Products”) against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product’s printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

**NEITHER ENVIROLOGIX NOR MONSANTO MAKE ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix’ sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix’ employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

**THIS WARRANTY IS EXCLUSIVE.** The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

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*Oster is a registered trademark of Sunbeam Products Inc.  
EnviroLogix, the EnviroLogix logo, and QuickStix are trademarks of EnviroLogix Inc.*

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#### License/Licencia/Licença

EnviroLogix has developed this kit using proprietary reagents as well as reagents licensed from Monsanto Company.

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#### For Technical Support Contact Us At:

**EnviroLogix**  
500 Riverside Industrial Parkway  
Portland, ME 04103-1486 USA

**Tel: (207) 797-0300**  
**Toll Free (U.S.): 866-408-4597**  
**Fax: (207) 797-7533**

e-mail: [info@envirologix.com](mailto:info@envirologix.com)

website: [www.envirologix.com](http://www.envirologix.com)

#### Para suporte técnico contactar:

**EnviroLogix do Brasil Diagnósticos Ltda.**  
Av. Emilio Marconato, 1000  
Galpão G24 Jd. Primavera  
CEP 13820-000, Jaguariúna, SP, Brasil

**Tel: +55 (19) 3307-8887**  
**Fax: +55 (19) 3307-8889**

[suportetecnico@envirologix.com](mailto:suportetecnico@envirologix.com)

E-mail: [vendas@envirologix.com](mailto:vendas@envirologix.com)

[www.envirologix.com.br](http://www.envirologix.com.br)

© EnviroLogix 2013

The test may be used qualitatively without the use of QuickScan. Presuming the low expressor data point, a visible pink test line will appear if the sample contains 0.3% Cry1Ac or more.

Alternatively, strips may be read visually at this point.

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