

### Set Contains:

- MB3 Extraction Buffer
- 1.5 mL clear micro-centrifuge tubes (50) for extraction
- 1.5 mL blue micro-centrifuge tubes (50) for dilution

### Materials Not Provided:

- Sterile Modified Buffered Peptone Water (mBPW), 2%
- Sterile collection sponge or swab for sampling surfaces. Collection sponge should be pre-moistened with Sterile Dey-Engley broth.
- Precision pipette(s) capable of delivering 25-1000 µL
- Pipette tips
- Incubator capable of 37±1°C
- Dry heat block capable of 95±1°C, with insert suitable for 1.5 mL tubes
- Vortexer
- Micro-centrifuge capable of 10,000 x g
- Timer
- Sterile Dey Engley (DE) culture broth (stainless steel surface testing only)
- Filter bags for culture

Catalog No. ACC-085

Part #11718

### Intended Use

- When used with DNABLE Kits for *Salmonella* (DF-026 and DF-126): This Set provides for extraction and detection of *Salmonella* DNA from a variety of matrices including dried pet food, wet pet food, stainless steel surfaces, poultry environmental swabs chicken carcass rinses, ground beef, powdered nonfat dry milk, lettuce, chocolate, and ground coriander. DNABLE *Salmonella* Supplement (Cat. No. XSALMD550) may be required for some matrices. Contact Technical Support for specific recommendations on other matrix types. This product is certified as a *Performance Tested Method*<sup>SM</sup> #041404 by the AOAC Research Institute as part of the DNABLE *Salmonella* assay for use for dry pet food, stainless steel surfaces and poultry environmental drag swabs.
- When used with DNABLE Kit for *Salmonella* Enteritidis (DF-056): This Set provides for extraction and detection of *Salmonella* Enteritidis DNA from a variety of matrices including boot swabs. DNABLE *Salmonella* Supplement (Cat. No. XSALMD500) may be required for some matrices. Contact Technical Support for specific recommendations on other matrix types.

### Intended User

DNABLE assays are designed to be simple and user friendly. They are designed for use by personnel with appropriate training in handling human pathogens, and in Microbiology and Molecular Assay techniques. Training specific to DNABLE assays will be provided by EnviroLogix; contact Technical Service or visit [envirologix.com/salmonella](http://envirologix.com/salmonella) for more information.

### How the Kit Works

An aliquot of MB3 buffer is added to a micro-centrifuge tube followed by a sample of enriched culture. The sample is concentrated by centrifugation and heated to enable lysis of *Salmonella* cells. A centrifugation and dilution step follows.

### Precautions and Notes

DNABLE is a highly sensitive assay. Therefore the following precautions are recommended to reduce the chance of sample contamination:

- Separated work areas are recommended for each of the following:
  - DNABLE culture
  - DNABLE sample preparation
  - DNABLE amplification and detection
- Clean the work station and pipettes with 10% bleach before and after use
- Do not reuse kit disposables
- Change pipette tips in between samples
- Wear gloves and change between handling of samples
- Avoid delays between sample preparation steps and between sample preparation and DNA amplification
- mBPW should be equilibrated at 37°C before use
- Filter bags should be used during enrichment to minimize particulates
- Enriched cultures should be mixed before sampling
- MB3 is stable for 1 year post manufacture when stored refrigerated at 4-8°C.

- **Safety:** *Salmonella* and *Salmonella* Enteritidis are pathogenic and classified as Biosafety Level 2 organisms. Personnel should be appropriately trained and should use personal protective equipment. Laboratories should follow appropriate local safety and environmental regulations and guidelines for containment and disposal as described in the Center for Disease Control and Prevention Manual, “Biosafety in Microbiological and Biomedical Laboratories” (link can be found here: [www.envirologix.com/useful-links](http://www.envirologix.com/useful-links)).

## Enrichment

### Notes

- Sample commodities to be tested according to your facility’s sampling plan.
- Composite samples according to procedures defined for your facility.
- Where required, prepare 2% supplemented media using DNable *Salmonella* Supplement (XSALMD500)
  - Prepare media according to manufacturer’s instructions and allow cooling to 30°C before adding supplement.
  - Using aseptic technique, add DNable *Salmonella* Supplement to media at a ratio of 2 to 100 and mix.
  - Once prepared, supplemented media may be stored refrigerated at 4-8°C for up to 9 months, protected from light.

Choose protocol based on sample type. Refer to the key listed with each sample type for applicability  
(**S** = *Salmonella* DF-026, **S+** = *Salmonella* Plus DF-126, **SE** = *Salmonella* Enteritidis DF-056):

<p><b>Dry pet food (S, S+)</b></p> <ol style="list-style-type: none"> <li>1. Suspend 25 g samples in 225 mL sterile mBPW.</li> <li>2. Suspend 375 g samples in either 1125 or 3375 mL mBPW. Use of filter culture bags is recommended.</li> <li>3. Allow the dry pet food to soften for 2-3 minutes, then stomach or otherwise homogenize the sample for 2 minutes (<math>\pm 30</math> seconds).</li> <li>4. Incubate 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>	<p><b>Wet pet food (S+)</b></p> <ol style="list-style-type: none"> <li>1. Suspend 375 g samples in 1125 mL mBPW. Use of filter culture bags is recommended.</li> <li>2. Allow the dry pet food to soften for 2-3 minutes, then stomach or otherwise homogenize the sample for 2 minutes (<math>\pm 30</math> seconds).</li> <li>3. Incubate 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>	<p><b>Poultry environmental drag swabs (S, S+, SE)</b></p> <ol style="list-style-type: none"> <li>1. Place Swab in 100 mL of mBPW.</li> <li>2. Shake the bag vigorously in an up-and-down motion at least ten times in an arc of approximately 30 cm over approximately 30 seconds</li> <li>3. Incubate 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>
<p><b>Stainless steel surface (S, S+)</b></p> <ol style="list-style-type: none"> <li>1. Pre-moisten collection sponge or swab with sterile DE broth.</li> <li>2. Swab the surface using horizontal and vertical motions.</li> <li>3. Place the sponge or swab in 75 mL of mBPW.</li> <li>4. Stomach by hand for 2 minutes (<math>\pm 30</math> seconds).</li> <li>5. Incubate suspended samples 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>	<p><b>Lettuce (S+)</b></p> <ol style="list-style-type: none"> <li>1. Suspend 375 g samples in 1125 mL mBPW. Use of filter culture bags is recommended.</li> <li>2. Stomach or otherwise homogenize the sample for 2 minutes (<math>\pm 30</math> seconds).</li> <li>3. Incubate 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>	<p><b>Chicken carcass (S, S+)</b></p> <ol style="list-style-type: none"> <li>1. Rinse whole bird as instructed in USDA-MLG (<a href="http://www.envirologix.com/useful-links">www.envirologix.com/useful-links</a>)</li> <li>2. Use 30 <math>\pm</math>0.5 mL of the sample rinse fluid obtained above for <i>Salmonella</i> analysis. Add 30 <math>\pm</math>0.5 mL of sterile mBPW and mix well.</li> <li>3. Incubate 20-22 hours at 37<math>\pm</math>1°C.</li> </ol>

(continued next page)

**Ground beef (S, S+)**

1. Suspend 375 g samples in 1125 modified tryptone soy broth, mTSB+n ([www.envirologix.com/useful-links](http://www.envirologix.com/useful-links)). Use of filter culture bags is recommended. Do not confuse mTSB+n media (MLG) with DNAble Modified TSB (dTSB) prepared using EnviroLogix *Salmonella* Supplement (XSALMD550).
2. Stomach or otherwise homogenize the sample for 2 minutes ( $\pm 30$  seconds).
3. Incubate 20-22 hours at  $37 \pm 1^\circ\text{C}$

*Contact Technical Service for additional matrix information and application guides for other matrices not listed herein*

**Sample Preparation**

1. Pre-heat a dry heat block to  $95^\circ\text{C}$ . Verify heat block is holding temperature with  $\pm 1.5^\circ$  using a simple thermometer.
2. Mix the culture before sampling. If testing dry pet food, squeeze the culture (375 g with 1125 mL of mBPW) to aid in sampling.
3. Transfer 1 mL of culture to a clear 1.5 mL micro-centrifuge tube supplied with the set.
4. Centrifuge the tube at  $10,000 \times g$  for 5 minutes ( $\pm 30$  seconds).
5. Remove the supernatant using caution to avoid disturbing the pellet. Leave a small volume remaining ( $\approx 100 \mu\text{L}$ ) if an obvious pellet is not observed.
6. Add  $100 \mu\text{L}$  of MB3 buffer and vortex to suspend the pellet.
7. Heat the sample in the heat block at  $95 \pm 1.5^\circ\text{C}$  for  $10 \pm 2$  minutes.
8. Centrifuge the tube at  $10,000 \times g$  for 5 minutes ( $\pm 30$  seconds).

*Choose next steps depending upon DNAble Assay:*

**S - Salmonella - DF-026**

9. Place **100  $\mu\text{L}$**  of MB3 into a blue 1.5 mL micro-centrifuge tube using a fresh pipette tip.
10. Transfer **25  $\mu\text{L}$**  of the supernatant from the second centrifugation into the MB3 and mix gently.

**S+ - Salmonella Plus - DF-126 and  
SE - Salmonella Enteritidis - DF-056**

9. Place **90  $\mu\text{L}$**  of MB3 into a blue 1.5 mL micro-centrifuge tube using a fresh pipette tip.
10. Transfer **10  $\mu\text{L}$**  of the supernatant from the second centrifugation into the MB3 and mix gently.

11. **5  $\mu\text{L}$**  of the sample from step 10 will be used in the subsequent DNAble reaction.

## LIMITED WARRANTY

EnviroLogix Inc. (“EnviroLogix”) warrants the products sold hereunder (“the Products”) against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product’s printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literatures are EnviroLogix’ sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix’ employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.



For Technical Support  
Contact Us At:

**EnviroLogix**  
500 Riverside Industrial Parkway  
Portland, ME 04103-1486 USA  
Tel: (207) 797-0300  
Toll Free: 866-408-4597  
Fax: (207) 797-7533

e-mail: [dnable@envirologix.com](mailto:dnable@envirologix.com)

Website: [www.envirologix.com](http://www.envirologix.com)

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

**THIS WARRANTY IS EXCLUSIVE.** The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

*EnviroLogix, the EnviroLogix logo, DNable, and the DNable logo are trademarks of EnviroLogix Inc.*

© EnviroLogix 2015